



STERIL TRANSPORT SWAB with CARY-BLAIR AGAR

INTENDED USE:

Transport swab with CARY BLAIR is used for collecting and storing microbe specimen. CARY BLAIR gel is effective in maintaining the viability of enteric pathogens, such as salmonella and shigella in rectal and fecal specimens.

COMPOSITION:

Transport swab with CARY BLAIR is composed of a swab, a transport tube, tube cap, label and CARY BLAIR gel.

CARY BLAIR gel Composition:

- Sodium thioglycollate
- Sodium hydrogen phosphate
- Calcium chloride
- Agar
- Sodium chloride

PRECAUTIONS:

- 1 For single use only, discards after use.
- 2 Store in dry and draughty warehouse.
- 3 Please do not use the product if the package is damaged.
4. In vitro diagnostic use only

STERILE CONDITION:

Sterilized by Gamma radiation sterility, the sterilization period is three years.

PACKAGING:

Packaging method: 1pc/blister, 1000pcs/carton.

Katalog Number: 08003

STORAGE CONDITIONS AND SHELF LIFE:

The product should be stored in dry and draughty warehouse. The period of validity is three years since manufacturing date, do not use the product beyond stated expiry date.